

Open modular ward trolleys - Configurations

Every ward trolley comes with a number of slots to fix the runners. Between every slot there is a gap of 1 inch. The chart below shows you how you can subdivide the trolleys.

Slots that are not in use are covered. To guarantee ease of use, we have left the 2 upper sections open.

These standard configurations are equipped with polyamide guides and a HPL worktop. Baskets and accessories are not included.

» Standard configurations

Configuration A Worktop height 33 1/2 inch 25 positions 25 1/2 x 17 3/4 x 34 inch	Configuration C Worktop height 39 1/2 inch 31 positions 25 1/2 x 17 3/4 x 40 inch	Configuration D Worktop height 39 1/2 inch 31 positions 13 3/4 x 17 3/4 x 40 inch	Configuration E Worktop height 42 1/2 inch 34 positions 25 1/2 x 17 3/4 x 43 inch

» Configurations on demand

Configuration B Worktop height 850 mm 25 positions 13 3/4 x 17 3/4 x 34 inch	Upon request The open modular nursing trolleys can also be closed with solid HPL walls. These nursing carts are available in different sizes.