

CSSD-SPD & Operating Theatre

CSSD-SPD & Operating Theatre

Sterisystem® Dry-Base® instrument trays	3
Steri-TOP lid	7
Steri-ID loop.....	11
Sterisystem® Dry-Base® Steri-ID instrument trays	13
CSSD Packing Tables.....	17
Stackers for Sterisystem® Perfo-Safe® wire baskets and shelves	23
Sterisystem® Perfo-Safe® wire baskets.....	25
Sterisystem® Perfo-Safe® wire shelves	29
Silicon mat	31
Sterisystem® open Stericart - 18006 series	33
Sterisystem® closed Stericart.....	37
Closed aluminum transport carts	43
Closed ISO stainless steel transport carts	45
UBeFlex® storage system	49
UBeTrack™ mobile storage system.....	59

Sterisystem® Dry-Base® instrument trays



Our Sterisystem® Dry-Base® instrument trays have been manufactured from highly polished AISI 304 grade stainless steel.

Our continuous quality control throughout the production process guarantees quality and dimensional stability. Therefore we offer a 10 year warranty on all our Sterisystem® Dry-Base® instrument trays. Thanks to their smooth and dimensionally stable base, you can rest assured that all your instruments will firmly stay in place. Careful attention has been paid to preventing perforations in packaging and great thought has been put into helping you optimise your sterilisation process.


Our Sterisystem® Dry-Base® instrument trays come with a perfectly finished top wire which reduces the chance of perforations to a minimum and a 10 year warranty.



»» Punch pattern

The fine, 5/64 inch holes made by the punch pattern guarantee a permeability of more than 50%, which can only benefit the washing process. As steam and water can readily penetrate during the sterilisation and washing process, an optimum result is guaranteed.

The sides protect the instruments against damage and have been designed to accommodate extra fine instruments.



»» Dry-Base®

The smooth, dry base completely traps any condensation moisture that may stick to the packaging. This ensures that condensation moisture can heat up or evaporate under the influence of radiation, convection or conduction.

Our Sterisystem® Dry-Base® instrument trays come with a perfectly finished, top wire which reduces the chance of perforations to a minimum.



»» Quality control tool

For quality-assurance purposes, all our trays are put through a special quality control tool. To make sure that the mesh is uniform in size, the base over every tray is meticulously checked, which is vital of course if we want to ensure that your instruments remain firmly in place.

Baskets that do not meet the dimensional or quality standards are rejected during quality control.



» Top wire

The upper wire is made from plasma that has been welded on to the sides of the tray. This results in a beautifully finished product that is devoid of any sharp edges and dispenses with the risk of perforations during unpacking. The upper wire also ensures that the top of the tray cannot become contaminated by dirt.



» Universal sizes

Both the DIN and ISO Sterisystem® Dry-Base® instrument trays come in universal sizes and fit perfectly into the Sterisystem® Perfo-Safe® wire trays.

Our Sterisystem® Dry-Base® instrument trays are fitted with 2 handles (except the DIN 1/4 and ISO 1/4 sizes) and have a load-bearing capacity of 33 lbs.



» Material

The Sterisystem® Dry-Base® instrument trays are manufactured from AISI 304 grade stainless steel and come with an electrolytically polished finish. The composition of the material is checked by an external knowledge centre.



» Shape

The Sterisystem® Dry-Base® instrument trays have been manufactured in such a way that you will never have to worry about holes or cavities that could harbour dirt.



» Bar code zone

For identification and traceability purposes, the Sterisystem® Dry-Base® instrument trays come with a bar code zone.

» Sterisystem® Dry-Base® instrument trays



» Type 60H - DIN Sterisystem® Dry-Base® instrument trays

Article Number	Description	GTIN	Ext. Dimensions
30006 00005	DIN 1/4	05404019100006	4 ^{23/32} x 9 ^{27/32} x 2 ^{33/64} inch
30006 00007	DIN 1/2	05404019100013	9 ^{29/64} x 9 ^{27/32} x 2 ^{33/64} inch
30006 00041	DIN 3/4	05404019100020	14 ^{11/64} x 9 ^{27/32} x 2 ^{33/64} inch
30006 00012	DIN 1/1	05404019100037	18 ^{57/64} x 9 ^{27/32} x 2 ^{33/64} inch



» Type 100H - DIN Sterisystem® Dry-Base® instrument trays

Article Number	Description	GTIN	Ext. Dimensions
30006 00020	DIN 1/2 H	05404019100044	9 ^{29/64} x 9 ^{27/32} x 3 ^{15/16} inch
30006 00011	DIN 1/1 H	05404019100051	18 ^{57/64} x 9 ^{27/32} x 3 ^{15/16} inch



» Type 60H - ISO Sterisystem® Dry-Base® instrument trays

Article Number	Description	GTIN	Ext. Dimensions
30006 00026	ISO 1/4	05404019100068	4 ^{23/32} x 11 ^{13/16} x 2 ^{33/64} inch
30006 00018	ISO 1/2	05404019100075	9 ^{29/64} x 11 ^{13/16} x 2 ^{33/64} inch
30006 00015	ISO 1/1	05404019100082	18 ^{57/64} x 11 ^{13/16} x 2 ^{33/64} inch



» Type 60H - ISO-2 Sterisystem® Dry-Base® instrument trays

Article Number	Description	GTIN	Ext. Dimensions
30006 00004	ISO-2 1/4	05404019100099	4 ^{17/32} x 12 ^{19/32} x 2 ^{33/64} inch
30006 00006	ISO-2 1/2	05404019100105	9 ^{1/16} x 12 ^{19/32} x 2 ^{33/64} inch
30006 00010	ISO-2 1/1	05404019100112	18 ^{7/64} x 12 ^{19/32} x 2 ^{33/64} inch



Steri-TOP lid



The Steri-TOP lid was designed with simple operation and guaranteed safety for use in a sterile environment in mind.

The instrument lids are made of highly polished AISI 304 grade stainless steel, in accordance with the same quality standards and requirements as our Sterisystem® Dry-Base® instrument trays.

The instrument tray is ergonomically opened by means of four silicone corners. In order to close the instrument tray, the four silicone corners are placed in the correct position and can easily be clicked onto the instrument tray.

The very small number of components, and the materials used, make the Steri-TOP lids the new standard for instrument tray lids, ensuring that no light instruments are blown out of the instrument tray during the sterilisation process.



» Steri-TOP lids sizes

All components are resistant to washing and autoclave sterilisation conditions (high temperatures, steam and vacuum) as well as to the most common chemicals used in hospitals.

The lids are available for DIN 1/4, 1/2, 3/4 and 1/1 and ISO 1/4, 1/2 and 1/1 baskets.



» Stackable

The ingenious design of the lids allows them to be stacked in storage when not in use. The silicone material of which the corners are made functions as an effective anti-slip layer.



» Matching instrument trays

The Steri-TOP lids perfectly match the Sterisystem® Dry-Base® instrument trays.

» Steri-TOP lid - DIN



» Steri-TOP lid DIN 1/4

Article number	Colour corners	GTIN	Ext. dimensions
30014 00021	Green	05404019100129	4 ^{23/32} x 9 ^{29/64} x 0 ^{43/64} inch
30014 00020	Yellow	05404019100136	4 ^{23/32} x 9 ^{29/64} x 0 ^{43/64} inch
30014 00022	Black	05404019100143	4 ^{23/32} x 9 ^{29/64} x 0 ^{43/64} inch
30014 00023	Blue	05404019100150	4 ^{23/32} x 9 ^{29/64} x 0 ^{43/64} inch
30014 00024	White	05404019100167	4 ^{23/32} x 9 ^{29/64} x 0 ^{43/64} inch
30014 00028	Red	05404019100174	4 ^{23/32} x 9 ^{29/64} x 0 ^{43/64} inch
30014 00056	Purple	05404019100181	4 ^{23/32} x 9 ^{29/64} x 0 ^{43/64} inch
30014 00057	Petrol blue	05404019100198	4 ^{23/32} x 9 ^{29/64} x 0 ^{43/64} inch



» Steri-TOP lid DIN 1/2

Article number	Colour corners	GTIN	Ext. dimensions
30014 00014	Green	05404019100204	9 ^{29/64} x 9 ^{27/32} x 0 ^{43/64} inch
30014 00013	Yellow	05404019100211	9 ^{29/64} x 9 ^{27/32} x 0 ^{43/64} inch
30014 00015	Black	05404019100228	9 ^{29/64} x 9 ^{27/32} x 0 ^{43/64} inch
30014 00016	Blue	05404019100235	9 ^{29/64} x 9 ^{27/32} x 0 ^{43/64} inch
30014 00017	White	05404019100242	9 ^{29/64} x 9 ^{27/32} x 0 ^{43/64} inch
30014 00027	Red	05404019100259	9 ^{29/64} x 9 ^{27/32} x 0 ^{43/64} inch
30014 00054	Purple	05404019100266	9 ^{29/64} x 9 ^{27/32} x 0 ^{43/64} inch
30014 00055	Petrol blue	05404019100273	9 ^{29/64} x 9 ^{27/32} x 0 ^{43/64} inch



» Steri-TOP lid DIN 3/4

Article number	Colour corners	GTIN	Ext. dimensions
30014 00066	Green	05404019102987	14 ^{11/64} x 9 ^{27/32} x 0 ^{43/64} inch

» Steri-TOP lid - ISO



» Steri-TOP lid DIN 1/1

Article number	Colour corners	GTIN	Ext. dimensions
30014 00009	Green	05404019100280	18 ^{57/64} x 9 ^{27/32} x 0 ^{43/64} inch
30014 00008	Yellow	05404019100297	18 ^{57/64} x 9 ^{27/32} x 0 ^{43/64} inch
30014 00010	Black	05404019100303	18 ^{57/64} x 9 ^{27/32} x 0 ^{43/64} inch
30014 00011	Blue	05404019100310	18 ^{57/64} x 9 ^{27/32} x 0 ^{43/64} inch
30014 00012	White	05404019100327	18 ^{57/64} x 9 ^{27/32} x 0 ^{43/64} inch
30014 00026	Red	05404019100334	18 ^{57/64} x 9 ^{27/32} x 0 ^{43/64} inch
30014 00052	Purple	05404019100341	18 ^{57/64} x 9 ^{27/32} x 0 ^{43/64} inch
30014 00053	Petrol blue	05404019100358	18 ^{57/64} x 9 ^{27/32} x 0 ^{43/64} inch



» Steri-TOP lid ISO 1/4

Article number	Colour corners	GTIN	Ext. dimensions
30014 00045	Green	05404019100365	4 ^{23/32} x 11 ^{13/16} x 0 ^{43/64} inch
30014 00044	Yellow	05404019100372	4 ^{23/32} x 11 ^{13/16} x 0 ^{43/64} inch
30014 00046	Black	05404019100389	4 ^{23/32} x 11 ^{13/16} x 0 ^{43/64} inch
30014 00047	Blue	05404019100396	4 ^{23/32} x 11 ^{13/16} x 0 ^{43/64} inch
30014 00048	White	05404019100402	4 ^{23/32} x 11 ^{13/16} x 0 ^{43/64} inch
30014 00049	Red	05404019100419	4 ^{23/32} x 11 ^{13/16} x 0 ^{43/64} inch
30014 00062	Purple	05404019100426	4 ^{23/32} x 11 ^{13/16} x 0 ^{43/64} inch
30014 00063	Petrol blue	05404019100433	4 ^{23/32} x 11 ^{13/16} x 0 ^{43/64} inch



» Steri-TOP lid ISO 1/2

Article number	Colour corners	GTIN	Ext. dimensions
30014 00038	Green	05404019100440	9 ^{29/64} x 11 ^{13/16} x 0 ^{43/64} inch
30014 00037	Yellow	05404019100457	9 ^{29/64} x 11 ^{13/16} x 0 ^{43/64} inch
30014 00039	Black	05404019100464	9 ^{29/64} x 11 ^{13/16} x 0 ^{43/64} inch
30014 00040	Blue	05404019100471	9 ^{29/64} x 11 ^{13/16} x 0 ^{43/64} inch
30014 00041	White	05404019100488	9 ^{29/64} x 11 ^{13/16} x 0 ^{43/64} inch
30014 00042	Red	05404019100495	9 ^{29/64} x 11 ^{13/16} x 0 ^{43/64} inch
30014 00060	Purple	05404019100501	9 ^{29/64} x 11 ^{13/16} x 0 ^{43/64} inch
30014 00061	Petrol blue	05404019100518	9 ^{29/64} x 11 ^{13/16} x 0 ^{43/64} inch

» Steri-TOP lid - DIN



» Steri-TOP lid ISO 1/1

Article number	Colour corners	GTIN	Ext. dimensions
30014 00031	Green	05404019100525	18 ^{57/64} x 11 ^{13/16} x 0 ^{43/64} inch
30014 00030	Yellow	05404019100532	18 ^{57/64} x 11 ^{13/16} x 0 ^{43/64} inch
30014 00032	Black	05404019100549	18 ^{57/64} x 11 ^{13/16} x 0 ^{43/64} inch
30014 00033	Blue	05404019100556	18 ^{57/64} x 11 ^{13/16} x 0 ^{43/64} inch
30014 00034	White	05404019100563	18 ^{57/64} x 11 ^{13/16} x 0 ^{43/64} inch
30014 00035	Red	05404019100570	18 ^{57/64} x 11 ^{13/16} x 0 ^{43/64} inch
30014 00058	Purple	05404019100587	18 ^{57/64} x 11 ^{13/16} x 0 ^{43/64} inch
30014 00059	Petrol blue	05404019100594	18 ^{57/64} x 11 ^{13/16} x 0 ^{43/64} inch

Steri-ID loop



The Steri-ID loop consists of a silicon loop in which an RFID tag is embedded.

The Steri-ID loop is a high quality complementary technology to the traditional process with barcode stickers for the tracing of the contents of sterilisation containers and Sterisystem® Dry-Base® instrument trays. Barcode stickers are easily damaged, become illegible and hence regularly need to be removed and replaced. The embedded tag however is resistant to washing and sterilization conditions in the autoclave (high temperatures, steam and vacuum) as well as to the most common chemicals used in hospitals. Therefore it can be used for over 400 cycles.

Because of their simple operation and sustainability, Steri-ID loops are a cost-efficient, time-saving and environmentally friendly alternative for the tracing of sterilization containers and instrument trays with barcode stickers.



» Steri-ID tag integrated

The built-in Steri-ID tag is based on ISO 15693 technology. The open platform allows you to choose from different suppliers of RFID writing and reading technologies.

The Steri-ID tag can be scanned through the sterilization packaging; hence no direct viewing is required to trace the contents during the entire sterilization process.



» Different colours

The Steri-ID loops are available in different colours. Herewith different standard surgical sets can be created and visually identified for example for gynecology, orthopedics,... or the sterilisation containers or Sterisystem® Dry-Base® instrument trays can be distinguished for the various campuses operating under the same CSSD.

» Steri-ID loop



» Steri-ID loop - HF

Article number	Colour
30012 00004	Green
30012 00003	Yellow
30012 00005	Black
30012 00006	Blue
30012 00007	White
30012 00028	Red
30012 00033	Purple
30012 00034	Petrol blue

Sterisystem® Dry-Base® Steri-ID instrument trays



Our Sterisystem® Dry-Base® instrument trays with integrated Steri-ID are made in accordance with the same quality criteria as our standard Sterisystem® Dry-Base® instrument trays.

The integrated Steri-ID HF RFID tag is a high quality complementary technology to the traditional process with barcode stickers used to keep track of the contents of Sterisystem® Dry-Base® instrument trays. These stickers are easily damaged, become illegible and as such regularly need to be removed and replaced. The embedded tag how-ever is resistant to washing and sterilization conditions in the autoclave (high temperatures, steam and vacuum) as well as to the most common chemicals used in hospitals. Therefore it can be used for over 400 cycles.

Because of their simple operation and sustainability, Sterisystem® Dry-Base® instrument trays with integrated Steri-ID tags are a cost-efficient, time-saving and environmentally friendly alternative to the traditional instrument trays with barcode stickers.



»» Different sizes and colours

The Steri-ID instrument tray is available in different standard sizes and colours. Herewith different standard surgical sets can be created and visually identified for example for gynecology, orthopedics,... or the surgical sets for the various campuses operating under the same CSSD can be distinguished.

The Sterisystem® Dry-Base® instrument trays are equipped with two handles (except the DIN 1/4 size) and have a load-bearing capacity of 33 lbs.



»» Integrated Steri-ID tag

The built-in Steri-ID tag is based on ISO 15693 technology. The open platform allows you to choose from different suppliers of RFID writing and reading technologies.

The Steri-ID tag can be scanned through the sterilization packaging; hence no direct viewing is required to trace the contents during the entire sterilization process.



»» Quality control

All our baskets are submitted to a quality control by a special control tool. The base of every basket is meticulously checked to make sure that the mesh is uniform in size. This is vital of course to ensure that the instruments remain firmly in place.


Baskets which do not comply with the dimensional or quality standards are rejected during the quality control.



» Punch pattern

The fine 5/64 inch holes of the punch pattern ensure a permeability of more than 50%, which is only beneficial for the washing process. As steam and water can readily penetrate during the washing and sterilization processes, a maximum result is guaranteed.

The sides protect the instruments against damage and are suitable for extra fine instruments.



» Dry-Base®

The smooth base remains in full contact with any condensation moisture left on the packaging. Subsequently radiation, convection or conduction can heat or evaporate this condensation moisture.



» Top wire

The top wire is plasma-welded on to the sides of the basket. This results in a perfectly finished product without any sharp edges, which reduces the risk of perforations during unpacking to a minimum.

The top wire is designed in such a way that no dirt can harbor between the sides and the top wire.



» Material and shape

The Sterisystem® Dry-Base® instrument trays are made off AISI 304 grade stainless steel and come with an electrolytically polished finish. The composition of the material is checked by an external knowledge center.


The Sterisystem® Dry-Base® instrument trays have been manufactured in such a way that no cavities are created which could harbor dirt.

» Sterisystem® Dry-Base® Steri-ID instrument trays



» Sterisystem® Dry-Base® Steri-ID instrument trays DIN 1/4

Article number	Tag colour	GTIN	Ext. dimensions
30012 00024	Green	05404019100600	4 ^{23/32} x 9 ^{27/32} x 2 ^{33/64} inch
Upon request	Yellow	05404019100617	4 ^{23/32} x 9 ^{27/32} x 2 ^{33/64} inch
Upon request	Black	05404019100624	4 ^{23/32} x 9 ^{27/32} x 2 ^{33/64} inch
Upon request	Blue	05404019100631	4 ^{23/32} x 9 ^{27/32} x 2 ^{33/64} inch
Upon request	White	05404019100648	4 ^{23/32} x 9 ^{27/32} x 2 ^{33/64} inch
Upon request	Red	05404019100655	4 ^{23/32} x 9 ^{27/32} x 2 ^{33/64} inch
Upon request	Purple	05404019100662	4 ^{23/32} x 9 ^{27/32} x 2 ^{33/64} inch
Upon request	Petrol blue	05404019100679	4 ^{23/32} x 9 ^{27/32} x 2 ^{33/64} inch



» Sterisystem® Dry-Base® Steri-ID instrument trays DIN 1/2

Article number	Tag colour	GTIN	Ext. dimensions
30012 00014	Green	05404019100686	9 ^{29/64} x 9 ^{27/32} x 2 ^{33/64} inch
Upon request	Yellow	05404019100693	9 ^{29/64} x 9 ^{27/32} x 2 ^{33/64} inch
Upon request	Black	05404019100709	9 ^{29/64} x 9 ^{27/32} x 2 ^{33/64} inch
Upon request	Blue	05404019100716	9 ^{29/64} x 9 ^{27/32} x 2 ^{33/64} inch
Upon request	White	05404019100723	9 ^{29/64} x 9 ^{27/32} x 2 ^{33/64} inch
Upon request	Red	05404019100730	9 ^{29/64} x 9 ^{27/32} x 2 ^{33/64} inch
Upon request	Purple	05404019100747	9 ^{29/64} x 9 ^{27/32} x 2 ^{33/64} inch
Upon request	Petrol blue	05404019100754	9 ^{29/64} x 9 ^{27/32} x 2 ^{33/64} inch



» Sterisystem® Dry-Base® Steri-ID instrument trays DIN 1/1

Article number	Tag colour	GTIN	Ext. dimensions
30012 00019	Green	05404019100761	18 ^{57/64} x 9 ^{27/32} x 2 ^{33/64} inch
Upon request	Yellow	05404019100778	18 ^{57/64} x 9 ^{27/32} x 2 ^{33/64} inch
Upon request	Black	05404019100785	18 ^{57/64} x 9 ^{27/32} x 2 ^{33/64} inch
Upon request	Blue	05404019100792	18 ^{57/64} x 9 ^{27/32} x 2 ^{33/64} inch
Upon request	White	05404019100808	18 ^{57/64} x 9 ^{27/32} x 2 ^{33/64} inch
Upon request	Red	05404019100815	18 ^{57/64} x 9 ^{27/32} x 2 ^{33/64} inch
Upon request	Purple	05404019100822	18 ^{57/64} x 9 ^{27/32} x 2 ^{33/64} inch
Upon request	Petrol blue	05404019100839	18 ^{57/64} x 9 ^{27/32} x 2 ^{33/64} inch

» Sterisystem® Dry-Base® Steri-ID instrument trays



» Replacement tags

Article number	Tag colour
30012 00009	Green
Upon request	Yellow
Upon request	Black
Upon request	Blue
Upon request	White
Upon request	Red
Upon request	Purple
Upon request	Petrol blue

CSSD Packing Tables



These height-adjustable packing tables have mainly been designed with people packing instrument sets in the central sterilisation unit in mind.

Our packing table allows your staff to work in an ergonomic fashion as the height can be adjusted to the type of task that needs to be performed and to people's anatomy and personal preferences. Thanks to the table's flexibility, staff can get on with their business in the most optimum working position.

The packing tables are shipped dismantled and should be assembled by your distributor's or your own technical staff.

The table's flexibility also means that workstations can be arranged to make the best possible use of the available space.



» Easily height adjustable

The height of the table can be set to anything between 30 ¹/₈ and 49 ⁵/₁₆ inches.

The height setting control comes with three different memory settings, allowing you to program a range of positions, which will prove very convenient if the same workstation is used by a number of people. The controls are located right underneath the worktop and are easy to reach.



» Flexibility

Aside from sockets, UTP and compressed air, the superstructure can also accommodate a range of accessories.

All CSSD Packing Tables come standard with 220V electrical supplies, if another type of electrical supply is required, this should be sourced locally.

Options and accessories marked with an * are intended for mounting during the production process, not afterwards.



» Worktop

The worktop is available in various sizes and finishes. Depending on the intended use, we can offer you a choice between stainless steel and HP laminate.

» CSSD Packing Tables - no upper structure fixture possible



» CSSD table with stainless steel worktop

Article Number	Dimensions worktop
17005 00086	70 ^{55/64} x 31 1/2 x 0 ^{25/32} inch
17005 00088	70 ^{55/64} x 39 3/8 x 0 ^{25/32} inch
17005 00090	78 ^{47/64} x 31 1/2 x 0 ^{25/32} inch
17005 00093	78 ^{47/64} x 39 3/8 x 0 ^{25/32} inch



» CSSD table with HP laminate worktop

Article Number	Dimensions worktop
17005 00192	70 ^{55/64} x 31 1/2 x 0 ^{33/64} inch
17005 00194	70 ^{55/64} x 39 3/8 x 0 ^{33/64} inch
17005 00196	78 ^{47/64} x 31 1/2 x 0 ^{33/64} inch
17005 00151	78 ^{47/64} x 39 3/8 x 0 ^{33/64} inch

» CSSD Packing Tables - suitable for upper structure components



» CSSD table stainless steel worktop for upper structure

Article Number	Dimensions worktop
17005 00082	70 ^{55/64} x 31 1/2 x 0 ^{25/32} inch
17005 00084	70 ^{55/64} x 39 3/8 x 0 ^{25/32} inch
17005 00046	78 ^{47/64} x 31 1/2 x 0 ^{25/32} inch
17005 00048	78 ^{47/64} x 39 3/8 x 0 ^{25/32} inch



» CSSD table with HP laminate worktop for upper structure

Article Number	Dimensions worktop
17005 00158	70 ^{55/64} x 31 1/2 x 0 ^{33/64} inch
17005 00165	70 ^{55/64} x 39 3/8 x 0 ^{33/64} inch
17005 00198	78 ^{47/64} x 31 1/2 x 0 ^{33/64} inch
17005 00149	78 ^{47/64} x 39 3/8 x 0 ^{33/64} inch

Photos not binding. The new generation of packing tables has drawer units with surface-mounted handles.

» Accessories for all types of CSSD Packing Tables



» Upper structure for baskets and trays *

Article Number	Description	Ext. Dimensions
17005 00022	With LED-lamp	55 ^{29/32} x 36 ^{7/32} inch
17005 00286	Without LED-lamp	55 ^{29/32} x 36 ^{7/32} inch



» Upper structure for cutting device *

Article Number	Description	Ext. Dimensions
17005 00116	With LED-lamp	55 ^{29/32} x 36 ^{7/32} inch
17005 00142	Without LED-lamp	55 ^{29/32} x 36 ^{7/32} inch
17005 00095	Cutting device	Length 47 ^{1/4} inch



» Shelf for stock of laminate sacs *

Article Number	Maximum weight	Dimensions shelf
17005 00131	22 lbs	15 ^{3/4} x 47 ^{63/64} x 0 ^{33/64} inch



» Sockets for upper structure *

Article Number	Description
17005 00023	3x socket
17005 00056	2x socket + switch (mounting always on the right-hand side)
17005 00057	2x socket + 2x UTP connection
17005 00058	1x Socket + 2x 2 UTP connections
17005 00160	1x socket + 2x UTP connection + switch (mounting always on the right-hand side)
17005 00148	Filling plate for electric sockets
17005 00190	Filling plate for electrical socket + switch (mounting always on the right-hand side)

Options and accessories marked with an * are intended for mounting during the production process, not afterwards.

» Upper structure components



» Stainless steel hook for the attachment of baskets

Article Number	Description
17005 00036	Oblique fixation of the wire basket
17005 00038	Straight fixation of the wire basket



» Sorting basket for upper structure

Article Number	Ext. Dimensions
15007 00043	5 ^{25/32} x 9 ^{1/16} x 5 ^{13/64} inch
15007 00042	8 ^{17/64} x 13 ^{25/32} x 5 ^{45/64} inch
15007 00044	8 ^{17/64} x 13 ^{25/32} x 7 ^{7/8} inch
17002 00023	Universal hook for the attachment of baskets to the upper structure



» Drawer unit including ABS modules with 3 drawers of 3^{15/16} inch high

Article Number	Description	Filling	Ext. Dimensions
17005 00250	Pre-mounted *	3 ABS modules	25 ^{43/64} x 18 ^{13/16} x 15 ^{63/64} inch
17005 00251	For mounting afterwards	3 ABS modules	25 ^{43/64} x 18 ^{13/16} x 15 ^{63/64} inch



» Flex arm with VESA 100 connection for attachment of a screen *

Article Number
14009 00019



» Inspection lamp with magnifying glass *

Article Number	Dimensions of the lens	Dioptre
30001 00034	Ø 4 ^{23/32} inch	5

Options and accessories marked with an * are intended for mounting during the production process, not afterwards.

» Upper structure components



» Inspection lamp with magnifying glass and adjustable light *

Article Number	Description	Dimensions of the lens	Dioptre
30001 00033	-	6 ^{3/8} x 4 ^{9/64} inch	3
On demand	With LED lighting	Ø 6 ^{1/32} inch	3,5



» Retractable drawer for Keyboard *

Article Number	Ext. Dimensions
17005 00167	18 ^{57/64} x 12 ^{1/64} inch




» Computerholder under table *

Article Number	Ext. Dimensions
17005 00070	7 ^{1/4} x 19 ^{31/64} x 19 ^{39/64} inch



» Wheels for limited displacement i.f.o. floor cleaning

Article Number
On request



» Compressed air connection

Article Number	Description
17005 00024	Fitted in 1 column of the upper structure

Options and accessories marked with an * are intended for mounting during the production process, not afterwards.



Stackers for Sterisystem®

Perfo-Safe® wire baskets and shelves



Made entirely of stainless steel, these spring-loaded trolleys are a real asset for your CSSD to help employees work more ergonomically.

Compared to storage on serving trolleys or in racks, you can gain a lot of space by storing Perfo-Safe® carriers in these stackers.

They are available for the ISO and DIN sizes of the Perfo-Safe® wire baskets and wire shelves.

Optionally, they can be equipped with a stainless steel push handle and wash tunnel-resistant wheels.

These high quality, stable and ergonomic stackers are a simple solution for the well-being and productivity of employees in the CSSD.



» Ergonomic functionality

Thanks to powerful springs tuned to the number and type of wire baskets or wire shelves in use, the baskets sink progressively downwards when loading. The same applies in the opposite direction when unloading. This allows CSSD employees to continuously work at ergonomic hip height (approx. 3.28 ft).



» Highly manoeuvrable and compact

The stackers have a reduced footprint and fit perfectly alongside Belintra's height-adjustable packing tables.

Thanks to the stainless steel 4 swivel castors, 2 of which are fitted with brakes at the front, the stackers are easy to move. Optionally, they can be equipped with wash tunnel-resistant castors.



» Robust and sturdy

The stackers have additional reinforcement of the base, legs and spring bracket to ensure long-term functionality. The number of springs and carrying capacity are matched to the type of wire baskets or wire grids that will be used in the trolley.

In addition, the stackers also have 4 bumpers in polypropylene.

» Stackers for Sterisystem® Perfo-Safe® wire baskets and shelves



» Stackers for Sterisystem® Perfo-Safe® wire baskets and shelves

Article Number	Description	Ext. Dimensions*
18019 00001	Stacker for DIN wire baskets	31 1/2 x 18 5/64 x 39 3/32 inch
18019 00002	Stacker for ISO wire baskets	31 1/2 x 22 1/64 x 39 3/32 inch
18019 xxxxxx	Stacker for DIN wire shelves	29 17/32 x 17 23/32 x 35 19/32 inch
18019 xxxxxx	Stacker for ISO wire shelves	29 17/32 x 21 21/32 x 35 19/32 inch
18019 00003	Stacker for DIN wire baskets with 4 wash tunnel resistant wheels	31 1/2 x 18 5/64 x 39 3/32 inch
18019 00004	Stacker for ISO wire baskets with 4 wash tunnel resistant wheels	31 1/2 x 22 1/64 x 39 3/32 inch

* Push handle excl.



» Stainless steel push handle for stackers

Article Number	Description
On request	Stainless steel push handle for DIN stackers
On request	Stainless steel push handle for ISO stackers

» Capacity

» Stacker for ISO wire baskets (18019 00002)

Article Number	Dimensions	Number
30005 00002 (Anna)	23 5/8 x 15 3/4 x 3 15/16 inch	28 baskets
30005 00004 (Alice)	23 5/8 x 15 3/4 x 5 1/8 inch	32 baskets
30005 00036 (Ava)	23 5/8 x 15 3/4 x 5 29/32 inch	30 baskets

» Stacker for DIN wire baskets (18019 00001)

Article Number	Dimensions	Number
30005 00028 (Abba)	23 5/8 x 11 13/16 x 4 21/64 inch	12 baskets
30005 00026 (Aline)	23 5/8 x 11 13/16 x 6 1/2 inch	8 baskets
30005 00034 (Amy)	23 5/8 x 11 13/16 x 1 31/32 inch	51 baskets

» Stacker for ISO wire shelves (18019 xxxxxx)

Article Number	Dimensions	Number
30003 00004 (Doris)	23 5/8 x 15 3/4 x 0 7/16 inch	70 shelves
30003 00007 (Delphine)	23 5/8 x 15 3/4 x 1 49/64 inch	58 shelves

» Stacker for DIN wire shelves (18019 xxxxxx)

Article Number	Dimensions	Number
30003 00003 (Deborah)	23 5/8 x 11 13/16 x 0 7/16 inch	70 shelves

Sterisystem® Perfo-Safe® wire baskets



Our Sterisystem® highly polished stainless steel Perfo-Safe® wire baskets, which are available in ISO (15 3/4 x 23 5/8 inch) and DIN (11 13/16 x 23 5/8 inch) sizes, are a safe and hygienic means to transport and store sterile packs.

As a sterilisation, transport and storage component, these stainless steel Sterisystem Perfo-Safe® wire baskets will prove to be an invaluable tool in your instrument set sterilisation process. As a matter of fact they reduce the number of touchpoints of the packaged instrument sets to just one after sterilisation. This reduces the risk for tears, cuts and holes in the packaging material.

These Sterisystem® Perfo-Safe® wire baskets are an indispensable part of sterile storage and transport.

Sterisystem® Perfo-Safe® wire baskets are invaluable if you want to ensure that your sterile instruments sets are transported and stored safely. Our Sterisystem® Perfo-Safe® wire baskets form the basis of the CSSD/SPD basket-in-basket system which reduces the risk for perforations.



» More sets per sterilisation cycle

With the medium and high DIN baskets you save quite some space in the autoclave. This thanks to the gain in width on the one hand, and the possibility of stacking them on top of each other on the other hand*.

The patented® design - with cutaways in the sides - guarantees that the instrument tray can be lifted ergonomically out of the wire basket.



» Compatible UBeFlex® DIN runners for storage

The DIN baskets are stored on top of the DIN runners in the notch provided for this purpose. This DIN runner is also suitable for the DIN Perfo-Safe® wire shelf.

This allows you to combine DIN wire baskets with DIN wire shelves according to your needs.



» Save space in transport

You can position two DIN baskets in the "W" (wide) direction after each other in the Belintra open or closed Stericarts. By doing so you can double the number of sets per level.



» Stackable* and nestable

Without the addition of extra handles, the baskets are not only stackable*, they are also nestable, making them more user-friendly and safer (no fear of fingers getting jammed).

By nesting the baskets which are not in use, you save costly space.



» Outstanding quality

The baskets have been manufactured in accordance with the most stringent quality standards. Baskets with welding burrs, gaps, wire distortions or discolourations in the welding work never see the light of day. Our baskets are devoid of any sharp edges or welding burrs which could cause tears in packaging.

For hygiene purposes, the upper edge features one single wire only instead of two.



» Increased patient safety

These Sterisystem Perfo-Safe® wire baskets are ideal if you wish to follow the basket-in-basket system to sterilise, transport and store your packaged Sterisystem Dry-Base® instrument trays. As the packaged instrument tray is stored in a wire basket, direct contact with the sterile instrument set is avoided.

» Sterisystem® Perfo-Safe® wire baskets



» Sterisystem® Perfo-Safe® Ergo-Stack wire baskets - DIN

Article Number	Name	GTIN	Ext. Dimensions
30005 00028	Abba	05404019100846	23 ^{5/8} x 11 ^{13/16} x 3 ^{15/16} inch
30005 00026	Aline	05404019100853	23 ^{5/8} x 11 ^{13/16} x 6 ^{1/2} inch



» Sterisystem® Perfo-Safe® wire basket - DIN

Article Number	Name	GTIN	Ext. Dimensions
30005 00001	Cindy*	05404019100907	23 ^{5/8} x 11 ^{13/16} x 1 ^{31/32} inch
30005 00034	Amy	05404019102291	23 ^{5/8} x 11 ^{13/16} x 1 ^{31/32} inch

* Only for use in UFlex 3 racks, not compatible with UBeFlex® racks.




» Sterisystem® Perfo-Safe® wire baskets - ISO

Article Number	Name	GTIN	Ext. Dimensions
30005 00002	Anna	05404019100860	23 ^{5/8} x 15 ^{3/4} x 3 ^{15/16} inch
30005 00004	Alice	05404019100877	23 ^{5/8} x 15 ^{3/4} x 5 ^{1/8} inch
30005 00036	Ava	05404019103069	23 ^{5/8} x 15 ^{3/4} x 5 ^{29/32} inch
30005 00018	Amanda	05404019100884	15 ^{3/4} x 11 ^{13/16} x 3 ^{15/16} inch
Upon request	Annetta	05404019100891	23 ^{5/8} x 15 ^{3/4} x 7 ^{7/8} inch



Sterisystem® Perfo-Safe® wire shelves



Our Sterisystem® highly polished Perfo-Safe® wire shelves are the perfect solution if you need to sterilise, transport and store diverse sterilisation packs. They are available in standard sizes : DIN (11^{13/16} x 23^{5/8} inch), ISO (15^{3/4} x 23^{5/8} inch) and for larger sterilisation packs in double DIN (23^{5/8} x 23^{5/8} inch).

Manufactured in accordance with the highest quality standards and finished to perfection, these wire trays meet the most stringent CSSD and OR requirements because the number of touchpoints of wrapped instrument trays is reduced to just one after sterilisation. This decreases the risk for tears, cuts and holes of the wrapping. These Sterisystem® Perfo-Safe® wire shelves are an indispensable part of sterile storage and transport.

Sterisystem® Perfo-Safe® wire shelves are an essential item if you wish to transport and store bulkier diverse sterilisation packs.



» Handles

The Sterisystem® Perfo-Safe® wire shelves have been fitted with 2 handles that have been seamlessly incorporated into the outer edge. These handles ensure that the wire trays are easy and ergonomic to manipulate.



» Welds

The Sterisystem® Perfo-Safe® wire shelves have been designed in such a way that any welds have been kept to an absolute minimum. They have been manufactured in accordance with the most stringent quality standards. Wire trays with welding burrs, gaps, wire distortions or discolourations in the welding work are not allowed.



» Rounded structure

Our Sterisystem® Perfo-Safe® wire shelves are devoid of any sharp edges or welding burrs that could cause tears, cuts and holes in packaging.

» Sterisystem® Perfo-Safe® wire shelves



» Sterisystem® Perfo-Safe® wire shelves


Article Number	Name	GTIN	Ext. Dimensions
30003 00003	Deborah	05404019100914	23 ^{5/8} x 11 ^{13/16} x 0 ^{7/16} inch
30003 00004	Doris	05404019100921	23 ^{5/8} x 15 ^{3/4} x 0 ^{7/16} inch
30003 00005	Dorothe	05404019100938	23 ^{5/8} x 23 ^{5/8} x 0 ^{7/16} inch
19031 00031*	XW	-	23 ^{5/8} x 37 ^{13/32} x 0 ^{7/16} inch



» Sterisystem® Perfo-Safe® wire shelf with raised edge

Article Number	Name	GTIN	Ext. Dimensions
30003 00007	Delphine	05404019100945	23 ^{5/8} x 15 ^{3/4} x 1 ^{49/64} inch
30003 00012	Demi	05404019103052	11 ^{13/16} x 15 ^{3/4} x 1 ^{49/64} inch

Silicon mat




The silicon mat is mainly used in the central sterilization department and is placed between the packaged instrument basket and the wire tray or wire basket.

Due to its excellent non-slip properties, the packaged instrument basket remains in place and does not move.


In combination with the Sterisystem® Perfo-Safe® wire shelves and baskets, the silicon mats are an indispensable part of sterile storage and transport.

The material of which the silicon mat is made is resistant to wear and to chemicals and bacteria, which makes it perfectly suitable for repeated use within the central sterilisation department.




» Anti-slip function

Since the packaged instrument basket does not move during transportation, any damage or perforation that could occur to the packaging material is avoided, guaranteeing sterility requirements within central sterilization department.



» Silicon


The material which the silicone mat is made from complies with requirements of the European Directive 2002/72/EC relating to plastic materials and articles intended to come into contact with foodstuffs, and is temperature resistant from -40 °F to 392 °F. This material is resistant to wear and to chemicals and bacteria, and corresponds to the highest degree of fire safety. Therefore, it is perfectly suitable for repeated use within the central sterilization department.



» Specific design

Its specific design ensures uniform heat distribution and perfect permeability of the steam within the autoclave. The standard colour is white. The thin version is appropriate for internal transportation, whereas the "heavy duty" 0^{15/64} inch thick version is suitable for external transportation.

» Silicon mat



» Silicon mat for Sterisystem® Perfo-Safe® wire baskets and wire shelves

Article Number	Description	Thickness	GTIN	Ext. Dimensions
30010 00004	Standard	0 ^{1/8} inch	05404019100952	21 ^{17/64} x 11 ^{1/32} inch
30010 00006	Heavy duty	0 ^{15/64} inch	05404019100969	23 ^{5/8} x 12 ^{3/4} inch

Sterisystem® open Stericart - 18006 series




The Sterisystem® open Stericart in stainless steel is designed for indoor transport of both instrument sets and consumables between the sterile processing department and the operating theatre.

This modular open transport cart comes both in an ISO size (23 ⁵/₈ x 15 ³/₄ inch) and in a double DIN size (23 ⁵/₈ x 23 ⁵/₈ inch). When using the double DIN version, the capacity per level is doubled if the establishment uses DIN transport baskets or trays.

The standard height of 59 ¹/₁₆ inch allows the user to see over the top of the cart during transport.

Sterisystem® open Stericarts are washing tunnel-proof, due to the high quality materials and the wheels which have been tested for resistance to washing tunnels.



» Hygiene

Underneath, the Sterisystem® open Stericarts are fitted with a base plate which slopes from the middle. This protects the transported goods against splashes of moisture or dirt.

This X-baseplate also ensures that no water remains behind during the washing and drying process.



» Washing tunnel-proof

The choice of high-quality, non-corrosive materials and components means that the Sterisystem® open Stericart can be cleaned in a washing tunnel.

Both the AISI 304 stainless steel body and the polyamide runners and wheels have been tested for resistance to washing tunnels.



» Unique handles

The characteristic yellow Belintra handles do not get as hot as stainless steel during washing, so the Sterisystem® open Stericart user can safely hold on to them at the end of the process.

The dual handles are installed vertically on either side so that they can be used ergonomically by both short and taller users.



» Noise-absorbing runners

The runners are made of plastic. Unlike stainless steel runners, the high quality plastic ensures smooth running and substantial noise reduction when sliding stainless steel transport baskets or trays or plastic modules in and out. The polyamide also absorbs noise during transport.

The runners have a built-in stop function at both ends. Every Sterisystem® open Stericart is fitted with 10 pairs of runners which results in 9 usable positions.



» User-friendly

Due to the functional height of 59 ^{1/16} inch, employees can still see over the top of the transport cart when moving forward. The high manoeuvrability of the wheels guarantees easy handling.

The open sides facilitate access on both sides, which certainly is an advantage in case of double storage per level.

» Versatile Modularity Tailored to Your Workflow

Tailor your cart to match your unique workflow and standards with our diverse range of open transport carts. Each model offers distinct modular features to optimize efficiency. Our modular open transport carts cater to the specific needs of your CSSD, providing adaptability without compromising efficiency.

Double DIN : Maximizing Capacity

The Double DIN (23 ^{5/8} x 23 ^{5/8} inch) model allows the insertion of two DIN transport baskets or shelves side by side on each level. This innovative design effectively doubles the capacity per level, accommodating up to 18 DIN sets. Additionally, the cart supports the use of 23 ^{5/8} x 23 ^{5/8} inch transport shelves, ISO transport baskets and shelves, as well as ISO PC and ABS modules for consumables.

ISO D-Type : Precision for ISO-Compatible Elements

Ideal for ISO-compatible elements, the ISO D-type model is optimized for ISO transport baskets and shelves, along with ISO PC and ABS modules for consumables. It excels in maintaining the integrity of ISO standards, though it does not support the insertion of DIN components.


ISO W-Type : A Blend of Flexibility

The ISO W-type model combines flexibility and adaptability. While it allows the insertion of DIN transport baskets or shelves (limited to one per level), it remains optimal for all ISO-compatible elements, including ISO transport baskets and shelves, as well as ISO PC and ABS modules for consumables.


2 x DIN : Doubling Efficiency

Like the Double DIN, the 2 x DIN version permits the insertion of two DIN transport baskets or shelves on each level in two separate columns, effectively doubling the capacity to accommodate up to 18 DIN sets. However, it's important to note that this model is not suitable for any ISO components.


» Possible filling per pair of runners (Double DIN)



2 x DIN wire shelves
23 ⁵/₈ x 11 ¹³/₁₆ inch
or
1 x ISO wire shelf
23 ⁵/₈ x 15 ³/₄ inch



2 x DIN wire baskets
23 ⁵/₈ x 11 ¹³/₁₆ inch
or
1 x ISO wire basket
23 ⁵/₈ x 15 ³/₄ inch



1 x double DIN wire shelf
23 ⁵/₈ x 23 ⁵/₈ inch



**1 x PC or ABS module
or basket**
23 ⁵/₈ x 15 ³/₄ inch

» Sterisystem® open Stericart - 18006 series



» Sterisystem® open Stericart - Double DIN

Article number	Runners	Ext. Dimensions
18006 00070	10 pair of fixed runners - 9 positions usable	31 ^{3/16} x 28 ^{5/32} x 61 ^{27/32} inch



» Sterisystem® open Stericart - ISO D-type

Article number	Runners	Ext. Dimensions
18006 00074	10 pair of fixed runners - 9 positions usable	23 ^{15/64} x 28 ^{5/32} x 61 ^{27/32} inch



» Sterisystem® open Stericart - ISO W-type

Article number	Runners	Ext. Dimensions
18006 00071	10 pair of fixed runners - 9 positions usable	31 ^{3/16} x 20 ^{9/32} x 61 ^{27/32} inch



» Sterisystem® open Stericart - 2 x DIN

Article number	Runners	Ext. Dimensions
18006 00092	2 x 10 pair of adjustable runners - 2 x 9 positions usable	32 ^{1/4} x 27 ^{3/32} x 65 ^{9/32} inch



» Sterisystem® open Stericart - 2 x DIN + shelf

Article number	Runners	Ext. Dimensions
18006 00097	2 x 7 pair of adjustable runners - 2 x 6 positions usable Open top compartment	32 ^{1/4} x 27 ^{3/32} x 65 ^{9/32} inch

Sterisystem® closed Stericart



Our Sterisystem® closed, washing-tunnel-resistant stainless steel Stericarts are manufactured based on the double DIN standard dimensions of 12 x 24 inch, with the idea of being able to use 2 DIN positions one after the other. The 2 inch extra depth at the back means that this transport cart can also be used for sterile transport of slightly longer sets. These Belintra Sterisystem® closed Stericarts are characterized by yellow runners, plastic bumpers, handlebars and latch.

The low model (43 inch) can double as a sterile field in the operating theatre. This case cart is a universal space-saving transport device. It is suitable for 2 DIN wire baskets or shelves or 1 ISO wire basket or shelf per level. The new "Extra-Wide" model can serve to transport packages up to 37 x 24 inch.

The Sterisystem® washing-tunnel-resistant and shock-resistant closed Stericarts with their unique space-saving design are perfect for safely transporting sterile instruments from and to the operating theatre and contaminated sets to the sterilizing unit.



» Wash tunnel resistant

The Sterisystem® closed case carts are designed in such a way to prevent any residual moisture after drying. The synthetic runners are also wash tunnel resistant*. Unlike stainless steel runners, the high quality plastic ensures smooth running and substantial noise reduction when sliding stainless steel transport baskets or trays or plastic modules in and out. Moreover, the synthetic runners protect wrapped instrument trays against perforations and absorb noise during transport.



» Ergonomic

The doors of these stainless steel transport carts swivel open to 270° and click onto the side walls with magnets when in the open position. This ensures that the cart is easy to access.

The high model still allows the user to see over the top of the cart during transport. The low model can also be used as a case cart with a sterile field in the operating theatre.



» Optimum loading and space-saving

The high case carts are equipped with 9 pairs of runners and the low case carts with 6 pairs. Thanks to the high number of positions, the loading is optimised without losing out on space.

Transporting the instruments on DIN Perfo-Safe® wire shelves or in DIN Perfo-Safe® wire baskets allows you to double the capacity per level.



» Safety

The cart can be sealed at the front with a breakable seal, to guarantee sterility from the sterilization unit to the operating theatre and also avoid contamination during the transport back.

It is recommended that the packed instrument baskets are placed on silicone anti-slip mats to prevent movement during transport and to reduce noise levels.



» Identification

The closed Stericarts feature a document holder on the side (A4) as standard to ensure that the carts and their contents can be identified from the outside.



» Castors

The transport carts are equipped with four double Tente® swivel castors with a diameter of 5 inch. Both wheels on the front have a brake and directional lock. The castors are equipped with a wheel hub on stainless steel bearings.

The carts are equipped with four bumpers to prevent damage to doors and walls.



» Cart washer resistant magnets

The Sterisystem® stainless steel transport cart is equipped with 4 sets (2 sets per door) of cart washer resistant neodymium magnets which are embedded in a silicone housing. Neodymium magnets are the strongest magnet on earth.

The design of the silicone housing allows for a perfect alignment of the magnets resulting in optimal connection power between the magnets and guaranteeing that the doors will remain open during the cleaning process in the cart washer.



» Handlebars

The low Sterisystem® transport carts are equipped with 1 horizontal yellow handlebar along each side, whereas the high carts are equipped with 2 vertical handlebars on each side to ensure smooth manoeuvrability.





» Status display

These transport trolleys are equipped with a status display where “green” indicates that the trolley contains clean, sterile goods. “Red” indicates that the cart is returning after use in the O.R.

Drying instructions: 8 minutes at 185°F in dry air. Dry air is air with a relative humidity of maximum 60% at an ambient temperature of 70°F. Converted this is a relative humidity of maximum 1,1% at 185°F.

» Sterisystem® closed Stericarts



» Sterisystem® closed Stericart - Low - Double DIN

Article number	Runners	Weight	Ext. Dimensions
19031 00039	1 column with 6 positions double DIN (23 ^{5/8} x 23 ^{5/8} inch)	128 lbs	28 ^{31/32} x 28 ^{31/32} x 43 ^{5/16} inch



» Sterisystem® closed Stericart - High - Double DIN

Article number	Runners	Weight	Ext. Dimensions
19031 00038	1 column with 9 positions double DIN (23 ^{5/8} x 23 ^{5/8} inch)	172 lbs	28 ^{31/32} x 28 ^{31/32} x 62 ^{63/64} inch



» Sterisystem® closed Stericart - XW


Article number	Columns	Weight	Ext. Dimensions
19031 00048	1 column with 6 positions for XW shelves (23 ^{5/8} x 37 ^{13/32} inch)	198 lbs	42 ^{43/64} x 28 ^{31/32} x 62 ^{63/64} inch



» Sterisystem® closed Stericart - 2 x DIN

Article number	Columns	Weight	Ext. Dimensions
19031 00044	2 columns with 6 DIN type positions (23 ^{5/8} x 11 ^{13/16} inch) 1 column with 2 positions for Dina shelf (25 ^{19/32} x 23 ^{5/8} inch) Art. nr. 19031 00047	181 lbs	30 ^{15/16} x 28 ^{31/32} x 62 ^{63/64} inch

» Carriers for use in the closed transport trolleys




2 x DIN wire shelves

24 x 12 inch
or

1 x ISO wire shelf

24 x 16 inch




2 x DIN wire baskets

24 x 12 inch
or


1 x ISO wire basket

24 x 16 inch




1 x double DIN wire shelf

24 x 24 inch



1 x PC or ABS module or basket


24 x 16 inch



XW wire shelf

24 x 37 inch

(Only for use in the Sterisystem®
transport cart - XW)



1 x Dina wire shelf

25^{19/32} x 24 inch



Closed aluminum transport carts



These modular trolleys in anodised aluminum are produced according to ISO and DIN measurements. The nylon guides, on which stainless steel as well as PC, ABS or polypropylene baskets or modules can rest, guarantee a smooth and noiseless transport. The distance between the guides is $4 \frac{31}{64}$ inch as standard.

Equipped with an extra bumper at the top of the trolley and using a lock, the aluminum transport trolley is a safe and ergonomic means of internal transport of sterile and non-sterile goods.

This transport cart is best cleaned manually. The use of aluminum trolleys in washing tunnels is not recommended, as the material will tarnish. For a washing tunnel resistant trolley, we recommend our range of stainless steel trolleys.

Our closed ISO aluminum transport carts allow you to safely and hygienically transport sterile and non-sterile goods or medication from your storage, sterilization unit or pharmacy to the relevant delivery points.



»» Tow bar mechanism

The base of the aluminum carts can be fitted with a tow bar to allow you to couple several carts together so that they can be operated as one single unit.

For enhanced manoeuvrability, both sides of the cart have been fitted with a push bar.



»» Removable ladders

If necessary, the ladders with the guides can be dismantled to allow for easy cleaning of the trolley.

The standard distance between the runners is $4 \frac{31}{64}$ inch.



»» Castors

The transport carts are delivered with $7 \frac{7}{8}$ inch diameter castors. These make for low rolling resistance and ensure that the transport cart is easy to manoeuvre.

Our 1-column, D-type transport carts have been fitted with smaller castors of $4 \frac{59}{64}$ inch in diameter.

» Closed aluminium transport carts - Overview

Basic version	19032 00001	19032 00002	19032 00003	19032 00004	19032 00005	19032 00006	19032 00007	19032 00008	19032 00009
All models inclusive of cylinder lock with key set, A4 document holder and brake on the swivel wheel									
Number of columns	ISO 1	ISO 1	ISO 2	ISO 2	ISO 3	ISO 3	DIN 2	ISO 3	ISO 2
Number of runners	12 pair	9 pair	2 x 12 pair	2 x 9 pair	3 x 12 pair	3 x 9 pair	2 x 9 pair	3 x 6 pair	2 x 6 pair
Usable positions (B = Baskets, S = Shelves)	12B 11S	9B 8S	2 x 12B 2 x 11S	2 x 9B 2 x 9S	3 x 12B 3 x 11S	3 x 9B 3 x 9S	2 x 9S	3 x 6B 3 x 6S	2 x 6B 2 x 6S
Document holder	✓	✓	✓	✓	✓	✓	✓	✓	✓
Drawbar	-	-	-	-	-	-	-	-	-
Model	High	Low	High	Low	High	Low	Low	Low	Low
Height in inches	65 ^{35/64}	59 ^{1/16}	68 ^{57/64}	59 ^{1/16}	68 ^{57/64}	59 ^{1/16}	59 ^{1/16}	59 ^{1/16}	59 ^{1/16}
Diameter wheels in inches	4 ^{59/64}	4 ^{59/64}	7 ^{7/8}	7 ^{7/8}	7 ^{7/8}	7 ^{7/8}	7 ^{7/8}	7 ^{7/8}	7 ^{7/8}
Load capacity	55 lbs per pair of runners - 330 lbs per column - 660 lbs per cart								
Filling									
ISO Perfo-Safe® wire shelves	✓	✓	✓	✓	✓	✓	-	✓	✓
ISO Perfo-Safe® wire baskets	✓	✓	✓	✓	✓	✓	-	✓	✓
ISO Perfo-Safe® wire shelf with raised edge (Delphine)	✓	✓	✓	✓	✓	✓	-	✓	✓
DIN Perfo-Safe® wire shelves	-	-	-	-	-	-	✓	-	-

Full option version	19032 00010	19032 00011	19032 00012	19032 00013	19032 00014	19032 00015	19032 00016	19032 00017	19032 00018
All models inclusive of cylinder lock with key set, A4 document holder and brake on the swivel wheel, rubber bumper at the top, tow bar and coupling hook									
Number of columns	ISO 1	ISO 1	ISO 2	ISO 2	ISO 3	ISO 3	DIN 2	ISO 3	ISO 2
Number of runners	12 pair	9 pair	2 x 12 pair	2 x 9 pair	3 x 12 pair	3 x 9 pair	2 x 9 pair	3 x 6 pair	2 x 6 pair
Usable positions (B = Baskets, S = Shelves)	12B 11S	9B 8S	2 x 12B 2 x 11S	2 x 9B 2 x 9S	3 x 12B 3 x 11S	3 x 9B 3 x 9S	2 x 9S	3 x 6B 3 x 6S	2 x 6B 2 x 6S
Document holder	✓	✓	✓	✓	✓	✓	✓	✓	✓
Drawbar	✓	✓	✓	✓	✓	✓	✓	✓	✓
Model	High	Low	High	Low	High	Low	Low	Low	Low
Height in inches	65 ^{35/64}	59 ^{1/16}	68 ^{57/64}	59 ^{1/16}	68 ^{57/64}	59 ^{1/16}	59 ^{1/16}	59 ^{1/16}	59 ^{1/16}
Diameter wheels in inches	4 ^{59/64}	4 ^{59/64}	7 ^{7/8}	7 ^{7/8}	7 ^{7/8}	7 ^{7/8}	7 ^{7/8}	7 ^{7/8}	7 ^{7/8}
Load capacity	55 lbs per pair of runners - 330 lbs per column - 660 lbs per cart								
Filling									
ISO Perfo-Safe® wire shelves	✓	✓	✓	✓	✓	✓	-	✓	✓
ISO Perfo-Safe® wire baskets	✓	✓	✓	✓	✓	✓	-	✓	✓
ISO Perfo-Safe® wire shelf with raised edge (Delphine)	✓	✓	✓	✓	✓	✓	-	✓	✓
DIN Perfo-Safe® wire shelves	-	-	-	-	-	-	✓	-	-

Closed ISO stainless steel transport carts



Our closed, washing-tunnel-resistant stainless steel transport carts are manufactured based on the ISO standard dimensions of 24 x 16 inch.

The carts are fitted with a lock with breakable seal, bumpers on the bottom, top and back, and identification holders on the front and sides.

This cart is ideal for transporting sterile sets from the external sterilization unit to the operating theatre. Used sets can also be transported from the hospital to the sterilization unit.

These ISO stainless steel washing-tunnel-resistant and shock-resistant transport carts are perfect for safely transporting sterile instruments from and to the operating theatre and contaminated sets to the sterilizing unit.



» Hygienic

The transport cart is washing-tunnel-resistant. The runners can be removed from the cart for optimum cleaning.

The base plate and top plate are sloped to prevent any residual moisture after drying.



» Door fixation

The doors of these stainless steel transport carts swivel open to 270° and are locked in the open position by a catch underneath the bumper. This ensures that the doors remain fixed in open position during the washing process.

The height (59 inch) allows the user to see over the top of the cart during transport.



» Optimal loading

The carts are equipped guides for each ISO column. These provide 9 useful positions depending on whether the instrument sets are being transported with ISO Perfo-Safe® wire shelves, wire baskets or wire shelves with raised edge (see comparative table).



» Safety

The carts can be sealed at the top with a breakable seal, to guarantee sterility from the sterilization unit to the operating theatre and also avoid contamination during the transport back.

It is recommended that the packed instrument baskets are placed on silicone anti-slip mats to prevent movement during external transport and to reduce noise levels.



» Identification

The transport carts feature a document holder on the front (A4) and side (A5) as standard to ensure that the carts and their contents can be identified from the outside.



» Castors

The transport carts are equipped with two swivel castors with a brake (left front and back) and two fixed castors (right front and back). The ISO 1 transport carts have a wheel diameter of 5 inch, while the ISO 2 transport carts have a wheel diameter of 6 ¹/₂ inch. Optionally, the cart can be equipped with directional fixation. This is a handy and ergonomic solution when logistics employees travel longer distances with **individual trolleys**. When connecting several carts in a train form, the combination of 2 fixed wheels and 2 swivel wheels is still a requirement.



» Drawbar

A number of standard configurations of the stainless steel transport carts are equipped with a drawbar that can be used to couple several transport carts together to form a 'train' (see comparative table).



» Status display

On request, transport trolleys can be equipped with a status display where "green" indicates that the trolley contains clean, sterile goods. "Red" indicates that the cart is returning after use in the O.R.

Model 19016 00115 is standard equipped with this status display.

» Closed ISO stainless steel transport carts - Overview

Article number	19016 00070	19016 00071	19016 00078	19016 00079
Number of columns	ISO 1	ISO 2	ISO 1	ISO 2
Number of runners	10 pair	2 x 10 pair	10 pair	2 x 10 pair
Usable positions	9	18	9	18
Document holder	✓	✓	✓	✓
Drawbar	-	-	✓	✓
Model	Low	Low	Low	Low
Height in inches	56	58	56	58
Diameter wheels in inches	5	6 1/2	5	6 1/2
Filling				
ISO Perfo-Safe® wire shelves	✓	✓	✓	✓
ISO Perfo-Safe® wire baskets	✓	✓	✓	✓
ISO Perfo-Safe® wire shelf with raised edge (Delphine)	-	-	-	-

Article number	19016 00076	19016 00077	19016 00080	19016 00081
Number of columns	ISO 1	ISO 2	ISO 1	ISO 2
Number of runners	9 pair	2 x 9 pair	9 pair	2 x 9 pair
Usable positions	9	18	9	18
Document holder	✓	✓	✓	✓
Drawbar	-	-	✓	✓
Model	Low	Low	Low	Low
Height in inches	56	58	56	58
Diameter wheels in inches	5	6 1/2	5	6 1/2
Filling				
ISO Perfo-Safe® wire shelves	✓	✓	✓	✓
ISO Perfo-Safe® wire baskets	-	-	-	-
ISO Perfo-Safe® wire shelf with raised edge (Delphine)	✓	✓	✓	✓

» Closed ISO stainless steel transport carts - Overview



Article number	19016 00115 - transport cart with shelf*
Number of shelves	1
Number of columns	ISO 2
Usable positions	7
Document holder	✓
Drawbar	-
Model	Low
Height in inches	59
Diameter wheels in inches	6 1/2
Filling	
ISO Perfo-Safe® wire shelves	7
ISO Perfo-Safe® wire baskets	7
ISO Perfo-Safe® wire shelf with raised edge (Delphine)	6-7

* Standard equipped with status display

UBeFlex® storage system



The UBeFlex® modular storage system was developed with very quick assembly, functionality, durability and stability in mind. The result is a unique storage system made from anodized aluminum and other high-quality, anti-static synthetic materials, which means it can be used in both sterile storage and non-sterile environments.

Thanks to the modular construction in DIN D-type, ISO D-type or ISO W-type columns, which may or may not be combined with shelves, it can be used in various departments of the healthcare institution, from CSA through pharmacy to kitchen and cold storage. Because of the unique design, it is an ergonomic and hygienic storage solution and can also be used as a “pass-through” rack.

With its revolutionary design, the new UBeFlex® storage system is very quickly assembled and is extremely flexible for optimum use of the available space.



» Easy mounting and cleaning

The columns, connection and connection sets, runners and shelves can be assembled quickly and easily without extra tools.

The UBeFlex® storage system can easily be cleaned and disinfected.



» Very modular with columns and shelving

UBeFlex® gives extensive possibilities for modular set-ups. With columns “Deep” or “Wide” 23 5/8 x 15 3/4 inch or DIN 23 5/8 x 11 13/16 inch in combination with shelves varying between 11 13/16 and 70 55/64 inch (11 13/16, 23 5/8, 35 7/16, 47 1/4, 59 1/16 or 70 55/64 inch) you can make optimum use of the available space.

For the shelving in D-type you can choose between full shelves in synthetic materials or open shelves in chromed steel. The open shelves can be fitted with dividers.



» Extremely strong

In column set-up, UBeFlex® has a load capacity of 55 lbs per pair of runners and in shelf set-up max. 33 lbs per shelf element. For details, see technical data sheet.



» Pass-through set-up

Because of the unique design, in which no reinforcement is required on the back, the UBeFlex® can also be used in a “pass-through” set-up and is thus accessible from both sides.

In case you want a back protection on the shelving, you can use additional shelving profiles to use as a closing profile to avoid packages sliding through.



» Stationary and movable versions

The UBeFlex® columns are available in a stationary version with 26 usable runner positions per column and a movable version with 24 usable runner positions per column. Unused positions can be covered with plastic caps.

The stationary version has height-adjustable feet to compensate for imperfections of the floor. The movable version is fitted with small wheels (Ø 1^{31/32} inch), so that it can be moved for floor cleaning.



» Optimum labeling

The shelves can be fitted with newly designed “click” label holders. These simply click onto any position on the front of the shelf and can be fitted with an insert. Because these label holders do not protrude above the edge of the shelf, you can take larger packets or containers from the shelves without a hitch.

There are specific label holders for the columns that can also simply be clicked onto the profiles. All versions can be found further back.



» Identification

As an option you can fit the UBeFlex® with coloured PVC strips (packs of 10) to replace the standard grey profiles to personalize your storage per department or, for example, to use colour coding for location allocation per discipline in the storage within your CSSD.

» New “Click” runners


Unique design of the plastic runners for easy assembly and quick movement if necessary. The runners come in packs of 10 pairs.



1. Position the runners in the required spot. Make sure that both pins have been inserted properly.




2. Push the runner against the upright and click the lip securely into place.



3. Make sure that the lip clicks properly over the edge of the upright.

» A set-up example


A UBeFlex® modular storage system module consists of 2 posts and a connection set. A connection set comprises of 4 connecting elements that interconnect the posts. To expand on this configuration you will need an extra post and connection set. Option wise, the sky is the limit.




» Combining columns of equal depth

Columns of equal depth can be combined ad infinitum. You can also interconnect columns of different widths. For the exact measurements, please refer to the next page.


» Example of a setup with runners



» Example of a setup with runners and shelves



» Example of a setup with shelves




In certain configurations, **spacer lockers** and **back reinforcement** are required. Be sure to ask your seller about this.

- The spacer lockers are included with the mobile posts, for the stationary posts these should be ordered separately if required.
- Back reinforcement is always required for configurations smaller than 3 columns and for some "heavy duty" applications.

» UBeFlex® Shelves

Every shelf of the UBeFlex® storage system consists of 2 profiles on which the boards are positioned. The entity forms the shelf on which goods can be placed.

Depending on the width of the shelf multiple boards and profiles of different lengths are required.



» Closing profiles

The profiles of the shelving can also be used as closing profiles at the back. Used as such the back side is blocked and packs can't fall off the shelving.



» Calculating the installation

By using the measures below you can calculate what the total width of your installation will be. You start with a **“Column”** which you then complete with one or more units **“Extension”**. ISO D-type and ISO DIN-type columns can be combined. Mixed installations with runners and shelves are also possible, as long as they have the same depth (D-type or W-type).

» Setup with runners

Type	Stationary	Movable (max. 5 columns coupled)
ISO D-type - Column	18 ^{31/32} inch	19 ^{17/32} inch
ISO D-type - Extension	17 ^{41/64} inch	17 ^{29/32} inch
DIN D-type - Column	15 ^{3/64} inch	15 ^{19/32} inch
DIN D-type - Extension	13 ^{45/64} inch	13 ^{31/32} inch
ISO W-type - Column	26 ^{27/32} inch	27 ^{13/32} inch
ISO W-type - Extension	25 ^{33/64} inch	25 ^{25/32} inch

» Setup with shelves

Type	Stationary	Movable (max. 5 columns coupled)
Shelving 11 ^{13/16} inch - Column	14 ^{3/8} inch	14 ^{59/64} inch
Shelving 11 ^{13/16} inch - Extension	13 ^{1/32} inch	13 ^{5/16} inch
Shelving 23 ^{5/8} inch - Column	26 ^{3/16} inch	26 ^{47/64} inch
Shelving 23 ^{5/8} inch - Extension	24 ^{27/32} inch	25 ^{1/8} inch
Shelving 35 ^{7/16} inch - Column	37 ^{63/64} inch	38 ^{35/64} inch
Shelving 35 ^{7/16} inch - Extension	36 ^{21/32} inch	36 ^{59/64} inch
Shelving 47 ^{1/4} inch - Column	49 ^{51/64} inch	50 ^{23/64} inch
Shelving 47 ^{1/4} inch - Extension	48 ^{15/32} inch	48 ^{47/64} inch
Shelving 59 ^{1/16} inch - Column	61 ^{39/64} inch	62 ^{11/64} inch
Shelving 59 ^{1/16} inch - Extension	60 ^{9/32} inch	60 ^{35/64} inch
Shelving 70 ^{55/64} inch* - Column	73 ^{27/64} inch	-
Shelving 70 ^{55/64} inch* - Extension	72 ^{3/32} inch	-

* The 70^{55/64} inch shelves cannot be used in the mobile version.



» UBeFlex® storage system



» Stationary posts

Article Number	Description	Depth	Height
16010 00091	ISO and DIN D-type	23 ^{5/8} inch	80 ^{5/16} inch
16010 00037	ISO W-type	15 ^{3/4} inch	80 ^{5/16} inch
16010 00066	Wall mounting (screw included)		
16010 00095	Floor mounting - stainless steel (including screws)		



» Posts with castors*

Article Number	Description	Depth	Height incl. castors
16010 00038	ISO and DIN D-type	23 ^{5/8} inch	77 ^{11/64} inch
16010 00039	ISO W-type	15 ^{3/4} inch	77 ^{11/64} inch

* Only for limited displacement in function of floor cleaning.

For open/closed internal and external transport, check www.belintra.com for the different solutions.



» Connection sets

Article Number	Description	Width
16010 00040	ISO D-type	15 ^{3/4} inch
16010 00041	ISO W-type	23 ^{5/8} inch
16010 00042	DIN D-type**	11 ^{13/16} inch
16010 00097	XS D-type	11 ^{1/32} inch



» Roof or bottom panel for columns with runners

Article Number	Description	Ext. Dimensions
16010 00030	ISO D-type	15 ^{3/4} x 23 ^{5/8} inch
16010 00031	ISO W-type	23 ^{5/8} x 15 ^{3/4} inch
16010 00065	DIN D-type	11 ^{13/16} x 23 ^{5/8} inch



» Cap for unused positions

Article Number	Description
99010 00020	100 pieces/bag

** For use in combination with D-type posts and DIN D-type runners. Filling with Sterisystem® Perfo-Safe® wire shelves and wire baskets DIN (11^{13/16} x 23^{5/8} inch).

» UBeFlex® storage system



» "Click"-runners (1 pair - left and right)

Article Number	Description	Type of runner	Length
16010 00011	ISO D-type	U+L runner	23 ^{5/8} inch
16010 00012	ISO W-type	U+L runner	15 ^{3/4} inch



» DIN D-type runners (1 pair - left and right)

Article Number	Description	Type of runner	Length
16010 00068	DIN D-type	L runner	23 ^{5/8} inch



» Oblique runners*

Article Number	Description	Fixation
14005 00005	ISO D-type oblique	Turn-pin
14005 00006	ISO W-type oblique	Turn-pin

* Upon request



» Catheter basket (without dividers)

Article Number	Description
16002 00020	ISO D-type - Equipped with telescopic runners and frame

» UBeFlex® storage system



» Catheter holder with adjustable hooks on melamine plate*

Article Number	Type	Description	Ext. Dimensions
1867580	ISO W-type	1 telescopic runner with 14 (2 x 7) hooks	23 ^{5/8} x 15 ^{3/4} inch
1867581	ISO W-type	2 telescopic runners with 14 (2 x 7) hooks	23 ^{5/8} x 15 ^{3/4} inch
1867582	ISO W-type	3 telescopic runners with 14 (2 x 7) hooks	23 ^{5/8} x 15 ^{3/4} inch
1867583	ISO D-type	1 telescopic runner with 20 (2 x 10) hooks	23 ^{5/8} x 15 ^{3/4} inch
1867584	ISO D-type	2 telescopic runners with 20 (2 x 10) hooks	23 ^{5/8} x 15 ^{3/4} inch

* Catheter holder only to be used in stationary racks.



» Catheter holder with fixed hooks on melamine plate*

Article Number	Type	Description	Ext. Dimensions
1867590	ISO W-type	1 telescopic runner with 14 (2 x 7) hooks	23 ^{5/8} x 15 ^{3/4} inch
1867591	ISO W-type	2 telescopic runners with 14 (2 x 7) hooks	23 ^{5/8} x 15 ^{3/4} inch
1867592	ISO W-type	3 telescopic runners with 14 (2 x 7) hooks	23 ^{5/8} x 15 ^{3/4} inch
1867593	ISO D-type	1 telescopic runner with 20 (2 x 10) hooks	23 ^{5/8} x 15 ^{3/4} inch
1867594	ISO D-type	2 telescopic runners with 20 (2 x 10) hooks	23 ^{5/8} x 15 ^{3/4} inch


* Catheter holder only to be used in stationary racks.



» UBeFlex® set of profiles for shelves

Article number	Description	
16010 00085	Set of profiles for shelf 11 ^{13/16} inch	Including hooks and caps, can also be used as a closing profile.
16010 00086	Set of profiles for shelf 23 ^{5/8} inch	
16010 00087	Set of profiles for shelf 35 ^{7/16} inch	
16010 00088	Set of profiles for shelf 47 ^{1/4} inch	
16010 00089	Set of profiles for shelf 59 ^{1/16} inch	
16010 00090	Set of profiles for shelf 70 ^{55/64} inch*	

* The 70^{55/64} inch shelves cannot be used in the mobile version.



» UBeFlex® shelves - excluding profiles

Article Number	Description	Ext. Dimensions
16010 00008	D-type full synthetic shelf	11 ^{13/16} inch W x 23 ^{5/8} inch D
16010 00009	W-type full synthetic shelf	11 ^{13/16} inch W x 15 ^{3/4} inch D
16010 00071	D-type open shelf in chromed steel	11 ^{13/16} inch W x 23 ^{5/8} inch D

» UBeFlex® storage system



» UBeFlex® D-type divider in chromed steel

Article Number

16010 00069



» Label holder for shelves - landscape

Article Number	Description	Ext. Dimensions
16010 00015	Label holder	-
15008 00021	Insert	2 ^{23/64} x 1 ^{37/64} inch



» Stepless label holder for adhesive labels for posts - portrait

Article Number	Description
16010 00121	Stepless label holder



» OptiScanBan® label holder for UBeFlex®/UBeTrack™ shelves

Article Number	Description
16010 00123	OptiScanBan® label holder



» Coloured PVC strips for posts

Article Number	Colour	Description	Length
16010 00076	Red	Conditioned per 10	79 ^{9/16} inch
16010 00077	Orange	Conditioned per 10	79 ^{9/16} inch
16010 00078	Blue	Conditioned per 10	79 ^{9/16} inch
16010 00079	Green	Conditioned per 10	79 ^{9/16} inch

» UBeFlex® storage system



» Telescopic runners with extension slide for containers* (MOQ: 100 pcs)

Article Number	Description	Max. weight	Ext. Dimensions
16002 00014	DIN D-type	33 lbs	11 ^{13/16} x 23 ^{5/8} inch
16010 00066	Wall fixation for posts		

* Only to be used when wall fixation is planned. To be mounted up to a maximum height of 39^{3/8} inch.



» Telescopic runners with extension slide for baskets*


Article Number	Description	Max. weight	Ext. Dimensions
16002 00008	ISO W-type	33 lbs	23 ^{5/8} x 15 ^{3/4} inch
16002 00012	ISO D-type	33 lbs	15 ^{3/4} x 23 ^{5/8} inch
16010 00066	Wall fixation for posts		

* Only to be used when wall fixation is planned. To be mounted up to a maximum height of 39^{3/8} inch.




» Spacer lockers

Article Number	Description
16010 00124	Included with the mobile posts, for the stationary posts to be ordered separately if required



» Spacer lockers in combination with roof/bottom panels

Article Number	Description
16010 00125	Included with the roof/bottom panels



» Back reinforcement kit

Article Number	Description
16010 00104	Required for certain configurations

UBeTrack™ mobile storage system



UBeTrack™ is a space-saving storage system with mobile racks.

The structure is very light so the mobile racks can be moved quickly. The result is a unique storage system made from anodized aluminium and high-quality, antistatic synthetic materials, which means it can be used in both sterile storage and non-sterile environments.

By positioning the runner system at the top good hygiene is guaranteed and it is easy to get past with transport trolleys.

The modular, space-saving UBeTrack™ storage system is a hygienic solution for small spaces with which you save up to 75% of corridor space and up to 33% of the total footprint.



» Easy assembly

The UBeTrack™ assembly is based on the UBeFlex® concept so it is easy and very quick to assemble. Since we use a single tracking system which is positioned on top of the unit and wheels on the bottom, there is no need for a track on the floor. Nothing to trip over!


The packages containing the parts are easy to carry and fit in a standard lift.



» Modular construction

ISO D-type and DIN D-type columns can be combined with each other and with D-type shelves of up to 47 1/4 inch. The same shelf-combinations must be fitted within one block however.

The maximum depth for working ergonomically is around 90 35/64 inch. For example, this corresponds with 5 ISO D-type columns or 6 DIN D-type columns whether combined with shelves or not.



» Hygienic

The tracking system, which is positioned at the top, has a protective cover over it preventing dust from falling on and accumulating in the track. This also allows for easy cleaning and disinfecting.

This top mounted track keeps floors clear and allows for easy cleaning of the floors and easy navigating through the unit with carts or trolleys.



» Design

With regard appearance the design fits in seamlessly with the existing UBeFlex® range so that both can be combined beautifully in one space.



» Castors

The wheels on the mobile shelves have very low rolling resistance.



» Extremely strong

In column set-up, UBeTrack™ has a load capacity of 55 lbs per pair of runners and in shelf set-up max. 30 lbs per shelf element. For details, see technical data sheet.



» Optimum labeling

The shelves can be fitted with newly designed “click” label holders. These simply click onto any position on the front of the shelf and can be fitted with an insert. Because these label holders do not protrude above the edge of the shelf, you can take larger packets or containers from the shelves without a hitch.

There are specific label holders for the columns that can also simply be clicked onto the profiles. All versions can be found further back.




» Identification


As an option you can fit the UBeTrack™ with coloured PVC strips (packs of 10) to replace the standard grey profiles to personalize your storage per department or, for example, to use colour coding for location allocation per discipline in the storage within your CSSD.

» Operation and construction of set-ups

The UBeTrack™ always consists of 1 full set-up that additional extension set-ups can be linked to.




One UBeTrack™ can consist of 2 or 3 mobile shelves. The total width can vary, see the illustration below. This can be increased in $0\frac{63}{64}$ inch increments up to a maximum of $9\frac{27}{32}$ inch which in turn alters the width of the walkway.



» UBeTrack™ mobile storage system

A UBeTrack™ installation consists of specific UBeTrack™ components completed with UBeFlex® parts.



» UBeFlex®/UBeTrack™ stationary post for outer racks

Article number	Description
16010 00091	ISO and DIN D-type



» UBeTrack™ mobile post for inner racks

Article number	Description
16011 00012	ISO and DIN D-type



» UBeTrack™ assembly kit

Article number	Description
16011 00015	1 per stationary rack



» UBeTrack™ track kit (2 to order per unit)

Article number	Description
16011 00013	59 ¹ / ₁₆ inch - For 2 mobile inner racks
16011 00014	72 ⁵³ / ₆₄ inch - For 3 mobile inner racks



» UBeTrack™ set of guide rollers for track kit

Article number	Description
16011 00018	1 per mobile rack

» UBeTrack™ mobile storage system



» UBeFlex® connection sets

Article Number	Description	Width
16010 00040	ISO D-type	15 ^{3/4} inch
16010 00042	DIN D-type	11 ^{13/16} inch



» UBeFlex® top or base panel for columns with runners

Article Number	Description	Ext. Dimensions
16010 00030	ISO D-type	15 ^{3/4} x 23 ^{5/8} inch
16010 00065	DIN D-type	11 ^{13/16} x 23 ^{5/8} inch

* Top or base panels are mandatory for columns with runners.



» UBeFlex® ISO D-type "Click"-runners (1 pair - left and right)

Article Number	Description	Type of runner	Length
16010 00011	ISO D-type	U+L runner	23 ^{5/8} inch



» UBeFlex® DIN D-type "Click"-runners (1 pair - left and right)

Article Number	Description	Type of runner	Length
16010 00068	DIN D-type	L runner	23 ^{5/8} inch



» UBeFlex® cap for unused positions

Article Number	Description
99010 00020	100 pieces/bag

» UBeTrack™ mobile storage system



» UBeFlex® D-type set for full synthetic shelves

Article Number	Width	Depth	Description
16010 00043	11 ¹³ / ₁₆ inch (1 x 11 ¹³ / ₁₆ inch)	23 ⁵ / ₈ inch	Including shelves, profiles, hooks and caps
16010 00044	23 ⁵ / ₈ inch (2 x 11 ¹³ / ₁₆ inch)	23 ⁵ / ₈ inch	



» UBeFlex® D-type shelves - excluding profiles

Article Number	Description	Ext. Dimensions
16010 00008	D-type full synthetic shelf	11 ¹³ / ₁₆ inch W x 23 ⁵ / ₈ inch D
16010 00071	D-type open shelf in chromed steel	11 ¹³ / ₁₆ inch W x 23 ⁵ / ₈ inch D



» UBeFlex® set of profiles for open shelves


Article number	Description	
16010 00085	Set of profiles for shelf 11 ¹³ / ₁₆ inch *	Including hooks and caps, can also be used as a closing profile.
16010 00086	Set of profiles for shelf 23 ⁵ / ₈ inch *	
16010 00087	Set of profiles for shelf 35 ⁷ / ₁₆ inch *	
16010 00088	Set of profiles for shelf 47 ¹ / ₄ inch *	



» UBeFlex® D-type divider in chromed steel

Article Number	Ext. Dimensions
16010 00069	11 ¹³ / ₁₆ x 23 ⁵ / ₈ inch

» UBeTrack™ mobile storage system



» Catheter holder with adjustable hooks on melamine plate*

Article Number	Type	Description	Ext. Dimensions
1867580	ISO W-type	1 telescopic runner with 14 (2 x 7) hooks	23 ^{5/8} x 15 ^{3/4} inch
1867581	ISO W-type	2 telescopic runners with 14 (2 x 7) hooks	23 ^{5/8} x 15 ^{3/4} inch
1867582	ISO W-type	3 telescopic runners with 14 (2 x 7) hooks	23 ^{5/8} x 15 ^{3/4} inch
1867583	ISO D-type	1 telescopic runner with 20 (2 x 10) hooks	23 ^{5/8} x 15 ^{3/4} inch
1867584	ISO D-type	2 telescopic runners with 20 (2 x 10) hooks	23 ^{5/8} x 15 ^{3/4} inch

* Catheter holder only to be used in stationary racks.



» Catheter holder with fixed hooks on melamine plate*

Article Number	Type	Description	Ext. Dimensions
1867590	ISO W-type	1 telescopic runner with 14 (2 x 7) hooks	23 ^{5/8} x 15 ^{3/4} inch
1867591	ISO W-type	2 telescopic runners with 14 (2 x 7) hooks	23 ^{5/8} x 15 ^{3/4} inch
1867592	ISO W-type	3 telescopic runners with 14 (2 x 7) hooks	23 ^{5/8} x 15 ^{3/4} inch
1867593	ISO D-type	1 telescopic runner with 20 (2 x 10) hooks	23 ^{5/8} x 15 ^{3/4} inch
1867594	ISO D-type	2 telescopic runners with 20 (2 x 10) hooks	23 ^{5/8} x 15 ^{3/4} inch

* Catheter holder only to be used in stationary racks.



» Label holder for shelves - landscape

Article Number	Description	Ext. Dimensions
16010 00015	Label holder	-
15008 00021	Insert	2 ^{23/64} x 1 ^{37/64} inch



» Stepless label holder for adhesive labels for posts - portrait

Article Number	Description
16010 00121	Stepless label holder

» UBeTrack™ mobile storage system



» OptiScanBan® label holder for UBeFlex®/UBeTrack™ shelves

Article Number	Description
16010 00123	OptiScanBan® label holder



» Coloured PVC strips for posts

Article Number	Colour	Description	Length
16010 00076	Red	Conditioned per 10	79 ⁹ / ₁₆ inch
16010 00077	Orange	Conditioned per 10	79 ⁹ / ₁₆ inch
16010 00078	Blue	Conditioned per 10	79 ⁹ / ₁₆ inch
16010 00079	Green	Conditioned per 10	79 ⁹ / ₁₆ inch



» Vertical handles for inner racks

Article Number	Description
<u>18002 00005</u>	Set of handles
16011 00017	UBeTrack™ assembly kit for vertical handles



» Bumper for UBeTrack™

Article Number	Description
16011 00019	Set of 4 synthetic bumpers incl. screws

Standard article (non-standard articles have a longer delivery time)



